

Floating Top Table

Suitable for use with our full range of x-ray generators, the floating top x-ray table makes accurate patient positioning quick and straightforward.

The radiolucent tabletop moves freely in all horizontal directions, allowing precise positioning with minimal patient disturbance. Once aligned, simply release the foot-operated lock to secure the tabletop quietly in position.

The cassette or DR panel tray is positioned beneath the tabletop, allowing a cassette to be changed without moving the patient. The x-ray generator remains automatically centred over the cassette or DR panel helping ensure consistent beam alignment for every exposure.

KEY FEATURES



Under surface self centering cassette and DR panel tray



Standard table top 1.6 m long



Optional 2.0 m table top



Suitable for all CR cassettes and DR panels



FLOATING TOP TABLE SPECIFICATION

Code	UN1007	Floating top table 1600 x 650mm
	UN1008	Floating top table 2000 x 650mm

